



Purpose-built to more effectively acquire tissue samples...

The Moray™ micro forceps is designed to enhance sampling from lesions that can occur within and outside of the gastrointestinal tract - leading to a more definitive diagnosis and targeted therapy.

The single-use Moray™ micro forceps features...

serrated jaws that effectively grab tissue

spring sheath allowing for use in tortuous positions

compatible with most 19g FNA needles

For more information on the Moray™ micro forceps, contact your local US Endoscopy representative today.

product number	description	sheath diameter (mm)	jaw diameter (mm)	jaw opening width (mm)	working length	units/box
BX00711889	Moray™ micro forceps	0.80	0.76	4.3	230cm	3

Distributed by:



listening...and delivering solutions®

5976 Heisley Road Mentor, OH 44060 USA

phone +1 440 / 639,4494 fax +1 440 / 639,4495

email global@usendoscopy.com